



LITERIDER 3-WHEEL

\$2,099.00

The LiteRider™ three-wheel scooter features up to a 12-mile range, a Delta Tiller standard, more attractive styling, and more vibrant colors. A modern tubular frame design supports its 300 lb. weight capacity and its completely wireless disassembly makes transporting a breeze.

SKU: GL111

Categories: [3-Wheel Scooters](#), [LiteRider Portable Mobility Scooters](#), [Mobility Scooters, Portable](#), [Travel Scooters](#)

ADDITIONAL INFORMATION

| | |
|---|-------------------------|
| Medicare Code | K0800 |
| Weight Capacity | 300 lbs |
| Drive Wheels | Rear |
| Maximum Speed | 5 mph |
| Operating Range | 12 miles |
| Ground Clearance | 3.75" |
| Ground Clearance to Transaxle/ Motor | 2.5" |
| Turning Radius | 39" |
| Batteries | 22 AH |
| Colors | Blue, Red |

| | |
|---|---|
| Features | 90 Amp Dynamic R Series Controller, Electro-Mechanical Brakes, Electronic Speed Control, Freewheel Mode |
| Length | 41.5" |
| Width | 22" |
| Height (ground to top of back of seat) | 37" - 39" |
| Ground to Top of Seat | 23.5" - 25.5" |
| Ground to Top of Deck | 6" |
| Deck to Top of Seat | 18" - 20" |
| Number Seat Height Positions | 3 |
| Size of Increments | 1" |
| Front Axle to Rear Axle | 30" |
| Weight w batteries | 125 lbs |
| Front Half Weight | 44 lbs |
| Rear Half Weight | 28 lbs |
| Seat weight w arms | 21 lbs |
| Batteries Weight | 29 lbs/pack |
| Heaviest Piece | 44 lbs |
| Front Tire | 9" |
| Rear tire | 9" |
| Rear Anti-Tippers | 2" |
| Standard Seat | Stadium Style |
| Seat Back Height - no headrest | 16" |
| Seat Width x Depth | 17" x 16" |
| Seat Color | Black Vinyl |
| Optional Seat | Stadium Style |
| Optional Seat Width x Depth | 20" x 18" |
| Optional Seat Color | Black Vinyl |

Warranty Drive Train 2 years

Warranty Batteries 1 year