

FDA Class II Medical Device



## LITERIDER 3-WHEEL

**\$2,099.00**

The LiteRider™ three-wheel scooter features up to a 12-mile range, a Delta Tiller standard, more attractive styling, and more vibrant colors. A modern tubular frame design supports its 300 lb. weight capacity and its completely wireless disassembly makes transporting a breeze.

---

**SKU:** GL111

**Categories:** [3-Wheel Scooters](#), [LiteRider](#)  
[Portable Mobility Scooters](#), [Mobility Scooters](#),  
[Portable](#), [Travel Scooters](#)

## ADDITIONAL INFORMATION

<b>Medicare Code</b>	K0800
<b>User Weight Capacity</b>	<a href="#">300 lbs</a>
<b>Drive Wheels</b>	Rear
<b>Maximum Speed</b>	5 mph
<b>Operating Range</b>	12 miles
<b>Ground Clearance</b>	3.75"
<b>Ground Clearance to Transaxle/ Motor</b>	2.5"
<b>Turning Radius</b>	39"
<b>Type Batteries</b>	22 AH
<b>Colors</b>	Blue, Red



<b>Features</b>	90 Amp Dynamic R Series Controller, Electro-Mechanical Brakes, Electronic Speed Control, Freewheel Mode
<b>Length</b>	41.5"
<b>Width</b>	22"
<b>Height (ground to top of back of seat)</b>	37" - 39"
<b>Ground to Top of Seat</b>	23.5" - 25.5"
<b>Ground to Top of Deck</b>	6"
<b>Deck to Top of Seat</b>	18" - 20"
<b>Number Seat Height Positions</b>	3
<b>Size of Increments</b>	1"
<b>Front Axle to Rear Axle</b>	30"
<b>Weight of Unit (assembled)</b>	125 lbs
<b>Front Half Weight</b>	44 lbs
<b>Rear Half Weight</b>	28 lbs
<b>Seat with Arms Weight</b>	21 lbs
<b>Batteries Weight</b>	29 lbs/pack
<b>Heaviest Piece</b>	44 lbs
<b>Front Tires</b>	9"
<b>Rear Tires</b>	9"
<b>Rear Anti-Tippers</b>	2"
<b>Standard Seat</b>	Stadium Style
<b>Seat Back Height - no headrest</b>	16"
<b>Seat Width x Depth</b>	18" x 16"
<b>Seat Color</b>	Black Vinyl
<b>Optional Seat</b>	Stadium Style
<b>Optional Seat Width x Depth</b>	20" x 18"
<b>Optional Seat Color</b>	Black Vinyl



**Warranty Drive Train** 2 years

**Warranty Batteries** 1 year