



LITERIDER 3-WHEEL

\$2,486.00

The LiteRider™ three-wheel scooter features up to a 12-mile range, a Delta Tiller standard, more attractive styling, and more vibrant colors. A modern tubular frame design supports its 300 lb. weight capacity and its completely wireless disassembly makes transporting a breeze.

SKU: GL111

Categories: [3-Wheel Scooters](#), [LiteRider Portable Mobility Scooters](#), [Mobility Scooters, Portable](#), [Travel Scooters](#)

ADDITIONAL INFORMATION

Medicare Code	K0800
Weight Capacity	300
Drive Wheels	Rear
Maximum Speed	5 MPH
Operating range	12 MILES
Ground Clearance	3.75"
Ground Clearance to Transaxle/ Motor	2.5"
Turning radius	39"
Batteries	22 AH
Colors	Blue, Red

Features	90 Amp Dynamic R Series Controller, Electro-Mechanical Brakes, Electronic Speed Control, Freewheel Mode
Length	41.5"
Width	22"
Height (ground to top of back of seat)	37" - 39"
Ground to Top of Seat	23.5" - 25.5"
Ground to Top of Deck	6"
Deck to Top of Seat	18" - 20"
Number Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	30"
Weight w batteries	125 LBS
Front Half Weight	44 LBS
Rear Half Weight	28 LBS
Seat weight w arms	21 LBS
Batteries Weight	29 lbs/pack
Heaviest Piece	44 LBS
Front Tire	9"
Rear tire	9"
Rear Anti-Tippers	2"
Standard Seat	Stadium Style
Seat Back Height - no headrest	16"
Seat Width x Depth	17" x 16"
Seat Color	Black Vinyl
Optional Seat	Stadium Style
Optional Seat Width x Depth	20" x 18"
Optional Seat Color	Black Vinyl



Warranty Drive Train 2 Year

Warranty Batteries 1 Year