



## LITERIDER 4-WHEEL

**\$2,299.00**

The four-wheel LiteRider™ scooter from Golden features 12.5-mile range, a Delta Tiller standard, more attractive styling and more vibrant colors.

---

**SKU:** GL141

**Categories:** [4-Wheel Scooters](#), [LiteRider Portable Mobility Scooters](#), [Mobility Scooters, Portable](#), [Travel Scooters](#)

## ADDITIONAL INFORMATION

<b>Medicare Code</b>	K0800
<b>Weight Capacity</b>	<a href="#">300 lbs</a>
<b>Drive Wheels</b>	Rear
<b>Maximum Speed</b>	5 mph
<b>Operating Range</b>	12.5 miles
<b>Ground Clearance Under Deck</b>	3.75"
<b>Ground Clearance to Transaxle/ Motor</b>	2.5"
<b>Turning Radius</b>	48"
<b>Batteries</b>	22 AH
<b>Colors</b>	Blue, Red
<b>Features</b>	90 Amp Dynamic R Series Controller, Electro-Mechanical Brakes, Electronic Speed Control, Freewheel Mode

<b>Charger</b>	Off Board
<b>Length</b>	43.25"
<b>Width</b>	22"
<b>Height (ground to top of back of seat)</b>	37" - 39"
<b>Ground to Top of Seat</b>	23.5" - 25.5"
<b>Ground to Top of Deck</b>	6"
<b>Deck to Top of Seat</b>	18" - 20"
<b>Number Seat Height Positions</b>	3
<b>Size of Increments</b>	1"
<b>Front Axle to Rear Axle</b>	30"
<b>Weight w batteries</b>	133 lbs
<b>Front Half Weight</b>	52 lbs
<b>Rear Half Weight</b>	28 lbs
<b>Seat weight w arms</b>	21 lbs
<b>Batteries Weight</b>	29 lbs/pack
<b>Heaviest Piece</b>	52 lbs
<b>Tires</b>	Solid/Low Profile
<b>Front Tire</b>	9"
<b>Rear tire</b>	9"
<b>Rear Anti-Tippers</b>	2"
<b>Standard Seat</b>	Stadium Style
<b>Seat Back Height - no headrest</b>	16"
<b>Seat Width x Depth</b>	17" x 16"
<b>Seat Color</b>	Black Vinyl
<b>Optional Seat</b>	Stadium Style
<b>Optional Seat Back Height - no headrest</b>	16"



<b>Optional Seat Width x Depth</b>	20" x 18"
<b>Optional Seat Color</b>	Black Vinyl
<b>Warranty Drive Train</b>	2 years
<b>Warranty Batteries</b>	1 year