

FDA Class II Medical Device



LITERIDER 4-WHEEL

\$2,299.00

The four-wheel LiteRider™ scooter from Golden features 12.5-mile range, a Delta Tiller standard, more attractive styling and more vibrant colors.

SKU: GL141

Categories: [4-Wheel Scooters](#), [LiteRider Portable Mobility Scooters](#), [Mobility Scooters, Portable](#), [Travel Scooters](#)

ADDITIONAL INFORMATION

Medicare Code	K0800
Weight Capacity	300 lbs
Drive Wheels	Rear
Maximum Speed	5 mph
Operating Range	12.5 miles
Ground Clearance Under Deck	3.75"
Ground Clearance to Transaxle/ Motor	2.5"
Turning Radius	48"
Batteries	22 AH
Colors	Blue, Red
Features	90 Amp Dynamic R Series Controller, Electro-Mechanical Brakes, Electronic Speed Control, Freewheel Mode

Charger	Off Board
Length	43.25"
Width	22"
Height (ground to top of back of seat)	37" - 39"
Ground to Top of Seat	23.5" - 25.5"
Ground to Top of Deck	6"
Deck to Top of Seat	18" - 20"
Number Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	30"
Weight w batteries	133 lbs
Front Half Weight	52 lbs
Rear Half Weight	28 lbs
Seat weight w arms	21 lbs
Batteries Weight	29 lbs/pack
Heaviest Piece	52 lbs
Tires	Solid/Low Profile
Front Tire	9"
Rear tire	9"
Rear Anti-Tippers	2"
Standard Seat	Stadium Style
Seat Back Height - no headrest	16"
Seat Width x Depth	17" x 16"
Seat Color	Black Vinyl
Optional Seat	Stadium Style
Optional Seat Back Height - no headrest	16"



Optional Seat Width x Depth	20" x 18"
Optional Seat Color	Black Vinyl
Warranty Drive Train	2 years
Warranty Batteries	1 year