

FDA Class II Medical Device

ALLY PRO

\$3,999.00



The Golden Ally Pro is a smart power wheelchair featuring a modern design with cutting-edge technologies. Folds easily & drives remotely.

FDA Class II Medical Devices designed by Golden are intended to assist persons with chronic disablements.

SKU: GP304

Categories: [Compact Power Chairs](#), [Foldable Power Wheelchairs](#), [Power Wheelchairs](#), [Travel Power Chairs](#)

ADDITIONAL INFORMATION

Medicare Code	Not Coded
User Weight Capacity	330 lbs
Drive Wheels	Rear
Maximum Speed	5 mph
Operating Range	13 miles
Ground Clearance	2.5"
Turning Radius	37"
Type Batteries	12 AH Lithium-Ion
Colors	Teal, Grey
Freewheel Mode	Yes



Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	40 Amp Remote
Length	39"
Width	25"
Height (ground to top of back of seat)	37"
Ground to Top of Seat	20.5"
Weight of Unit (assembled)	75 lbs.
Weight of Battery	5 lbs.
Tires (standard)	Pneumatic
Rear Tires	10"
Casters	10" Omni-Directional
Anti-Tip Wheels	2"
Standard Seat	Stadium Style
Back Height (no headrest)	19"
Standard Seat Width x Depth	17"x16"
Seat Color	Black
Footrest Height Adjustable	Yes
Footrest Angle Adjustable	Yes
Warranty Frame	5 Years
Warranty Drive Train	13 Months
Warranty Electronics	13 Months
Warranty Batteries	12 Months