



ALLY PRO

\$3,749.00

The Golden Ally Pro is a smart power wheelchair featuring a modern design with cutting-edge technologies. Folds easily & drives remotely.

SKU: GP304

Categories: [Compact Power Chairs](#), [Favorites](#), [Foldable Power Wheelchairs](#), [Power Wheelchairs](#), [Travel Power Chairs](#)

ADDITIONAL INFORMATION

Class II Medical Device	Yes
Medicare Code	Not Coded
Weight Capacity	330 lbs
Drive Wheels	Rear
Maximum Speed	5 mph
Operating Range	13 miles
Ground Clearance	2.5"
Turning Radius	37"
Batteries	12 AH Lithium-Ion
Colors	Teal, Grey
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes

Charger	Off Board
Controller	40 Amp Remote
Length	39"
Width	25"
Width Between Arms	20" Max
Height (ground to top of back of seat)	37"
Ground to Top of Seat	20.5"
Weight w batteries	75 lbs.
Weight of Battery	5 lbs.
Heaviest Piece	69.6 lbs.
Drive Wheels	10" Pneumatic
Casters	10" Omni-Directional
Anti-Tip Wheels	2"
Standard Seat Width x Depth	17x16"
Footrest Height Adjustable	Yes
Increments	1"
Warranty Frame	5 Years
Warranty Drive Train	13 Months
Warranty Electronics	13 Months
Warranty Batteries	12 Months