

FDA Class II Medical Device



## ALLY SPORT

**\$4,249.00**

The Golden Ally Sport is a smart, luxurious power wheelchair featuring a modern design with cutting-edge technologies.

FDA Class II Medical Devices designed by Golden are intended to assist persons with chronic disabilities.

---

**SKU:** GP305

**Categories:** [Compact Power Chairs](#), [Foldable Power Wheelchairs](#), [Power Wheelchairs](#), [Travel Power Chairs](#)

## ADDITIONAL INFORMATION

<b>Medicare Code</b>	Not Coded
<b>User Weight Capacity</b>	<a href="#">265 lbs</a>
<b>Drive Wheels</b>	Rear
<b>Maximum Speed</b>	4 mph
<b>Operating Range</b>	9.5 miles
<b>Ground Clearance</b>	2.5"
<b>Turning Radius</b>	30"
<b>Type Batteries</b>	10 AH Lithium-Ion
<b>Colors</b>	White, Black Carbon Fiber Pattern
<b>Freewheel Mode</b>	Yes
<b>Electronic Speed Control</b>	Yes



<b>Electro-Mechanical Brakes</b>	Yes
<b>Charger</b>	Off Board
<b>Controller</b>	40 Amp Remote
<b>Length</b>	42"
<b>Width</b>	25"
<b>Height (ground to top of back of seat)</b>	37"
<b>Ground to Top of Seat</b>	22.5"
<b>Weight of Unit (assembled)</b>	92 lbs.
<b>Weight of Battery</b>	5 lbs.
<b>Tires (standard)</b>	Pneumatic
<b>Casters</b>	10" Omni-Directional
<b>Anti-Tip Wheels</b>	2"
<b>Standard Seat</b>	Stadium Style
<b>Back Height (no headrest)</b>	18"
<b>Standard Seat Width x Depth</b>	17" x 17"
<b>Seat Color</b>	Black, Grey
<b>Footrest Height Adjustable</b>	No
<b>Footrest Angle Adjustable</b>	Yes
<b>Warranty Frame</b>	5 Years
<b>Warranty Drive Train</b>	13 Months
<b>Warranty Electronics</b>	13 Months
<b>Warranty Batteries</b>	12 Months