

FDA Class II Medical Device



ALLY SPORT

\$4,249.00

The Golden Ally Sport is a smart, luxurious power wheelchair featuring a modern design with cutting-edge technologies.

FDA Class II Medical Devices designed by Golden are intended to assist persons with chronic disablements.

SKU: GP305

Categories: [Compact Power Chairs](#), [Favorites](#), [Foldable Power Wheelchairs](#), [Power Wheelchairs](#), [Travel Power Chairs](#)

Tags: [New](#)

ADDITIONAL INFORMATION

Medicare Code	Not Coded
User Weight Capacity	265 lbs
Drive Wheels	Rear
Maximum Speed	4 mph
Operating Range	9.5 miles
Ground Clearance	2.5"
Turning Radius	30"
Type Batteries	10 AH Lithium-Ion
Colors	White, Black Carbon Fiber Pattern
Freewheel Mode	Yes
Electronic Speed Control	Yes



Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	40 Amp Remote
Length	42"
Width	25"
Height (ground to top of back of seat)	37"
Ground to Top of Seat	22.5"
Weight of Unit (assembled)	92 lbs.
Weight of Battery	5 lbs.
Tires (standard)	Pneumatic
Casters	10" Omni-Directional
Anti-Tip Wheels	2"
Standard Seat	Stadium Style
Back Height (no headrest)	18"
Standard Seat Width x Depth	17" x 17"
Seat Color	Black, Grey
Footrest Height Adjustable	No
Footrest Angle Adjustable	Yes
Warranty Frame	5 Years
Warranty Drive Train	13 Months
Warranty Electronics	13 Months
Warranty Batteries	12 Months