



## ALLY

**\$3,629.00**

Experience freedom with the new Golden Ally™, our new power chair designed to provide the strength and solid feel of a traditional chair yet remaining light weight and portable enough to be the perfect travel partner.

FDA Class II Medical Devices designed by Golden are intended to assist persons with chronic disablements.

---

**SKU:** GP303

**Categories:** [Compact Power Chairs](#), [Foldable Power Wheelchairs](#), [Power Wheelchairs](#), [Travel Power Chairs](#)

## ADDITIONAL INFORMATION

<b>Medicare Code</b>	Not Coded
<b>User Weight Capacity</b>	<a href="#">330 lbs</a>
<b>Drive Wheels</b>	Rear
<b>Maximum Speed</b>	5 mph
<b>Operating Range</b>	9 miles
<b>Ground Clearance</b>	2.5"
<b>Turning Radius</b>	28.5"
<b>Type Batteries</b>	12AH 24V Lithium-Ion
<b>Colors</b>	Teal, White, Black



<b>Freewheel Mode</b>	Yes
<b>Electronic Speed Control</b>	Yes
<b>Electro-Mechanical Brakes</b>	Yes
<b>Charger</b>	Off Board
<b>Controller</b>	40 Amp Remote
<b>Length</b>	41"
<b>Width</b>	24.5"
<b>Height (ground to top of back of seat)</b>	36"
<b>Ground to Top of Seat</b>	21.5"
<b>Weight of Unit (assembled)</b>	56 lbs
<b>Weight of Battery</b>	5 lbs.
<b>Tires (standard)</b>	Pneumatic
<b>Rear Tires</b>	10"
<b>Casters</b>	6"
<b>Rear Anti-Tippers</b>	2"
<b>Standard Seat</b>	Stadium Style
<b>Back Height (no headrest)</b>	18"
<b>Seat Width x Depth</b>	17" x 16"
<b>Seat Color</b>	Black
<b>Footrest Height Adjustable</b>	Yes
<b>Footrest Angle Adjustable</b>	Yes
<b>Warranty Frame</b>	5 Years
<b>Warranty Drive Train</b>	13 months
<b>Warranty Electronics</b>	13 months
<b>Warranty Batteries</b>	1 year