



ALLY

\$3,629.00

Experience freedom with the new Golden Ally™, our new power chair designed to provide the strength and solid feel of a traditional chair yet remaining light weight and portable enough to be the perfect travel partner.

SKU: GP303

Categories: [Compact Power Chairs](#), [Foldable Power Wheelchairs](#), [Power Wheelchairs](#), [Travel Power Chairs](#)

ADDITIONAL INFORMATION

Class II Medical Device	Yes
Medicare Code	Not Coded
Weight Capacity	330 lbs
Drive Wheels	Rear
Maximum Speed	5 mph
Operating Range	9 miles
Ground Clearance	2.5"
Turning Radius	28.5"
Batteries	12 AH Lithium-Ion
Colors	Green, White, Black
Freewheel Mode	Yes



Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	40 Amp Remote
Length	41"
Width	24.5"
Height (ground to top of back of seat)	36"
Ground to Top of Seat	21.5"
Weight w batteries	56 lbs
Weight of Battery	5 lbs.
Drive Wheels	10" Pneumatic
Standard Seat Width x Depth	17" x 16"
Footrest Height Adjustable	Yes
Increments	1"
Warranty Frame	5 Years
Warranty Drive Train	13 Months
Warranty Electronics	13 Months
Warranty Batteries	12 Months