

**ALLY****\$3,629.00**

Experience freedom with the new Golden Ally™, our new power chair designed to provide the strength and solid feel of a traditional chair yet remaining light weight and portable enough to be the perfect travel partner.

---

**SKU:** GP303

**Categories:** [Compact Power Chairs](#), [Favorites](#), [Foldable Power Wheelchairs](#), [Power Wheelchairs](#), [Travel Power Chairs](#)

**ADDITIONAL INFORMATION**

<b>Class II Medical Device</b>	Yes
<b>Medicare Code</b>	Not Coded
<b>Weight Capacity</b>	<a href="#">330 lbs</a>
<b>Drive Wheels</b>	Rear
<b>Maximum Speed</b>	5 mph
<b>Operating Range</b>	9 miles
<b>Ground Clearance</b>	2.5"
<b>Turning Radius</b>	28.5"
<b>Batteries</b>	12 AH Lithium-Ion
<b>Colors</b>	Green, White, Black
<b>Freewheel Mode</b>	Yes



<b>Electronic Speed Control</b>	Yes
<b>Electro-Mechanical Brakes</b>	Yes
<b>Charger</b>	Off Board
<b>Controller</b>	40 Amp Remote
<b>Length</b>	41"
<b>Width</b>	24.5"
<b>Height (ground to top of back of seat)</b>	36"
<b>Ground to Top of Seat</b>	21.5"
<b>Weight w batteries</b>	56 lbs
<b>Weight of Battery</b>	5 lbs.
<b>Drive Wheels</b>	10" Pneumatic
<b>Standard Seat Width x Depth</b>	17" x 16"
<b>Footrest Height Adjustable</b>	Yes
<b>Increments</b>	1"
<b>Warranty Frame</b>	5 Years
<b>Warranty Drive Train</b>	13 Months
<b>Warranty Electronics</b>	13 Months
<b>Warranty Batteries</b>	12 Months