

FDA Class II Medical Device



ALLY

\$3,629.00

Experience freedom with the new Golden Ally™, our new power chair designed to provide the strength and solid feel of a traditional chair yet remaining light weight and portable enough to be the perfect travel partner.

FDA Class II Medical Devices designed by Golden are intended to assist persons with chronic disablements.

SKU: GP303

Categories: [Compact Power Chairs](#), [Foldable Power Wheelchairs](#), [Power Wheelchairs](#), [Travel Power Chairs](#)

ADDITIONAL INFORMATION

Medicare Code	Not Coded
User Weight Capacity	330 lbs
Drive Wheels	Rear
Maximum Speed	5 mph
Operating Range	9 miles
Ground Clearance	2.5"
Turning Radius	28.5"
Type Batteries	12AH 24V Lithium-Ion
Colors	Teal, White, Black



Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	40 Amp Remote
Length	41"
Width	24.5"
Height (ground to top of back of seat)	36"
Ground to Top of Seat	21.5"
Weight of Unit (assembled)	56 lbs
Weight of Battery	5 lbs.
Tires (standard)	Pneumatic
Rear Tires	10"
Casters	6"
Rear Anti-Tippers	2"
Standard Seat	Stadium Style
Back Height (no headrest)	18"
Seat Width x Depth	17" x 16"
Seat Color	Black
Footrest Height Adjustable	Yes
Footrest Angle Adjustable	Yes
Warranty Frame	5 Years
Warranty Drive Train	13 months
Warranty Electronics	13 months
Warranty Batteries	1 year