



## BUZZABOUT

**\$2,399.00**

Enjoy full-time luxury on the go in the new BuzzAbout! Perfect for travel, the BuzzAbout power wheelchair is in a class all its own in our popular BuzzAround series of mobility products. The BuzzAbout has a stylish captain's seat with contoured support and higher back for longer seating tolerance and comfort.

FDA Class II Medical Devices designed by Golden are intended to assist persons with chronic disablements.

---

**SKU:** GP164

**Categories:** [Compact Power Chairs](#), [Power Wheelchairs](#), [Travel Power Chairs](#)

**Tags:** [New](#)

## ADDITIONAL INFORMATION

<b>Medicare Code</b>	Not Coded
<b>User Weight Capacity</b>	<a href="#">300 lbs</a>
<b>Drive Wheels</b>	Rear
<b>Maximum Speed</b>	4 mph
<b>Operating Range</b>	11 miles
<b>Ground Clearance</b>	2.5"
<b>Turning Radius</b>	28.5"
<b>Type Batteries</b>	22 AH



<b>Standard Colors</b>	Red, Blue
<b>Freewheel Mode</b>	Yes
<b>Electronic Speed Control</b>	Yes
<b>Electro-Mechanical Brakes</b>	Yes
<b>Charger</b>	Off Board
<b>Controller</b>	40 Amp
<b>Length</b>	36"
<b>Width</b>	23"
<b>Height (ground to top of back of seat)</b>	37.5" - 39.5"
<b>Height (ground to top of back of headrest)</b>	36.75" - 38.75"
<b>Ground to Top of Seat</b>	21.5" - 23.5"
<b>Number Seat Height Positions</b>	3
<b>Size of Increments</b>	1"
<b>Weight of Unit (assembled)</b>	124 lbs
<b>Front Half Weight</b>	29 lbs
<b>Rear Half Weight</b>	33 lbs
<b>Seat with Arms Weight</b>	34 lbs
<b>Batteries Weight</b>	28 lbs.
<b>Tires (standard)</b>	Flat Free / Non Marking Solid
<b>Rear Tires</b>	9"
<b>Casters</b>	6"
<b>Rear Anti-Tippers</b>	2"
<b>Standard Seat</b>	Highback
<b>Back Height (no headrest)</b>	18"
<b>Back Height with Headrest</b>	30.5" - 34.5"
<b>Seat Width x Depth</b>	17" x 16"
<b>Seat Color</b>	Black Vinyl



<b>Footrest Height Adjustable</b>	No
<b>Footrest Angle Adjustable</b>	Yes
<b>Warranty Frame</b>	5 Years
<b>Warranty Drive Train</b>	13 months
<b>Warranty Electronics</b>	13 months
<b>Warranty Batteries</b>	1 year