

**BUZZABOUT****\$2,399.00**

Enjoy full-time luxury on the go in the new BuzzAbout! Perfect for travel, the BuzzAbout power wheelchair is in a class all its own in our popular BuzzAround series of mobility products. The BuzzAbout has a stylish captain's seat with contoured support and higher back for longer seating tolerance and comfort.

---

**SKU:** GP164**Categories:** [Compact Power Chairs](#), [Power Wheelchairs](#), [Travel Power Chairs](#)**Tags:** [New](#)**ADDITIONAL INFORMATION**

<b>Length</b>	36"
<b>Width</b>	23"
<b>Front Tire</b>	N/A
<b>Rear tire</b>	9"
<b>Center Tire</b>	N/A
<b>Casters</b>	6"
<b>Weight Capacity</b>	<a href="#">300 lbs</a>
<b>Maximum Speed</b>	4 mph
<b>Ground Clearance</b>	2.5"
<b>Operating Range</b>	11 miles



<b>Turning Radius</b>	28.5"
<b>Batteries</b>	22 AH
<b>Weight w batteries</b>	124 lbs
<b>Front Half Weight</b>	29 lbs
<b>Rear Half Weight</b>	33 lbs
<b>Seat weight w arms</b>	34 lbs
<b>Heaviest Piece</b>	34 lbs
<b>Base Weight</b>	N/A
<b>Warranty on Electrical Parts</b>	13 months
<b>Features</b>	Adjustable Armrests, Adjustable Headrest, Adjustable Seat Posts, Center Wheel Drive, Electronic Speed Control, Footrest Extension, Freewheel Mode
<b>Warranty Batteries</b>	1 year