

FDA Class II Medical Device



BUZZABOUT

\$2,399.00

This compact, portable chair features motion-sensing LED lights on the footplate and rear handle that activate with joystick movement for added safety and visibility.

The BuzzAbout also offers convenient offboard charging, easy disassembly for transport, and excellent maneuverability—indoors or out.

SKU: GP164

Categories: [Compact Power Chairs](#), [Power Wheelchairs](#), [Travel Power Chairs](#)

ADDITIONAL INFORMATION

Medicare Code	Not Coded
User Weight Capacity	300 lbs
Drive Wheels	Rear
Maximum Speed	4 mph
Operating Range	11 miles
Ground Clearance	2.5"
Turning Radius	28.5"
Type Batteries	22 AH
Standard Colors	Red, Blue
Freewheel Mode	Yes



Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	40 Amp
Length	36"
Width	23"
Height (ground to top of back of seat)	37.5" - 39.5"
Height (ground to top of back of headrest)	36.75" - 38.75"
Ground to Top of Seat	21.5" - 23.5"
Number Seat Height Positions	3
Size of Increments	1"
Weight of Unit (assembled)	124 lbs
Front Half Weight	29 lbs
Rear Half Weight	33 lbs
Seat with Arms Weight	34 lbs
Batteries Weight	28 lbs.
Tires (standard)	Flat Free / Non Marking Solid
Rear Tires	9"
Casters	6"
Rear Anti-Tippers	2"
Standard Seat	Highback
Back Height (no headrest)	18"
Back Height with Headrest	30.5" - 34.5"
Seat Width x Depth	17" x 16"
Seat Color	Black Vinyl
Footrest Height Adjustable	No
Footrest Angle Adjustable	Yes



Warranty Frame	5 Years
Warranty Drive Train	13 months
Warranty Electronics	13 months
Warranty Batteries	1 year