

**BUZZABOUT****\$2,399.00**

Enjoy full-time luxury on the go in the new BuzzAbout! Perfect for travel, the BuzzAbout power wheelchair is in a class all its own in our popular BuzzAround series of mobility products. The BuzzAbout has a stylish captain's seat with contoured support and higher back for longer seating tolerance and comfort.

SKU: GP164**Categories:** [Compact Power Chairs](#), [Power Wheelchairs](#), [Travel Power Chairs](#)**Tags:** [New](#)**ADDITIONAL INFORMATION**

Length	36"
Width	23"
Front Tire	N/A
Rear tire	9"
Center Tire	N/A
Casters	6"
Weight Capacity	300 lbs
Maximum Speed	4 mph
Ground Clearance	2.5"
Operating Range	11 miles

Turning Radius	28.5"
Batteries	22 AH
Weight w batteries	124 lbs
Front Half Weight	29 lbs
Rear Half Weight	33 lbs
Seat weight w arms	34 lbs
Heaviest Piece	34 lbs
Base Weight	N/A
Warranty on Electrical Parts	13 months
Features	Adjustable Armrests, Adjustable Headrest, Adjustable Seat Posts, Center Wheel Drive, Electronic Speed Control, Footrest Extension, Freewheel Mode
Warranty Batteries	1 year