

FDA Class II Medical Device



BUZZAROUND CARRYON

\$3,199.00

The Buzzaround CarryOn power chair is super portable with a 21" turning radius and center wheel drive for optimal control and maneuverability.

SKU: GP130

Categories: [Center Wheel Power Chairs](#), [Compact Power Chairs](#), [Foldable Power Wheelchairs](#), [Power Wheelchairs](#), [Travel Power Chairs](#)

ADDITIONAL INFORMATION

Medicare Code	N/A
User Weight Capacity	300 lbs
Drive Wheels	Center
Maximum Speed	3.85 mph
Operating Range	12.5 miles
Ground Clearance	2"
Turning Radius	21"
Type Batteries	12 AH Lithium-Ion
Colors	Red, Blue
Optional Color Shroud Kits	Grey, White, Yellow
Freewheel Mode	Yes
Electro-Mechanical Brakes	Yes



Charger	Off Board
Controller	40 Amp
Length	38"
Width	24"
Height (ground to top of back of seat)	34.25"
Ground to Top of Seat	21"
Weight of Unit Assembled	74 lbs.
Weight of Base	46 lbs.
Weight of Seat	22 lbs.
Weight of Battery	6 lbs.
Drive Wheels	7.75"
Standard Seat	Stadium Style
Back Height - No Headrest	19.3"
Seat Width x Depth	17" x 16"
Color	Black
Warranty Frame	5 Years
Warranty Drive Train	13 Months
Warranty Electronics	13 Months
Warranty Batteries	12 Months