



BUZZAROUND CARRYON

\$3,199.00

The Buzzaround CarryOn power chair is super portable with a 21" turning radius and center wheel drive for optimal control and maneuverability.

SKU: GP130

Categories: [Center Wheel Power Chairs](#), [Compact Power Chairs](#), [Favorites](#), [Foldable Power Wheelchairs](#), [Power Wheelchairs](#), [Travel Power Chairs](#)

Tags: [New](#)

ADDITIONAL INFORMATION

Medicare Code	N/A
Weight Capacity	300 lbs
Drive Wheels	Center
Maximum Speed	3.85 mph
Operating Range	12.5 miles
Ground Clearance	2"
Turning Radius	21"
Batteries	12 AH Lithium-Ion
Colors	Red, Blue
Freewheel Mode	Yes
Electro-Mechanical Brakes	Yes
Charger	2A, Off Board

Controller	40 Amp
Length	38"
Width	24"
Height (ground to top of back of seat)	34.25"
Ground to Top of Seat	21"
Weight of Unit Assembled	74 lbs.
Weight of Base	46 lbs.
Weight of Seat	22 lbs.
Weight of Battery	6 lbs.
Drive Wheels	7.75"
Standard Seat Width x Depth	18" x 17"
Back Height - No Headrest	15.5"
Color	Black
Warranty Frame	5 Years
Warranty Drive Train	13 Months
Warranty Electronics	13 Months
Warranty Batteries	12 Months