



## BUZZAROUND CARRYON

**\$3,199.00**

The Buzzaround CarryOn power chair is super portable with a 21" turning radius and center wheel drive for optimal control and maneuverability.

---

**SKU:** GP130

**Categories:** [Center Wheel Power Chairs](#), [Compact Power Chairs](#), [Foldable Power Wheelchairs](#), [Power Wheelchairs](#), [Travel Power Chairs](#)

**Tags:** [New](#)

## ADDITIONAL INFORMATION

<b>Medicare Code</b>	N/A
<b>Weight Capacity</b>	<a href="#">300 lbs</a>
<b>Drive Wheels</b>	Center
<b>Maximum Speed</b>	3.85 mph
<b>Operating Range</b>	12.5 miles
<b>Ground Clearance</b>	2"
<b>Turning Radius</b>	21"
<b>Batteries</b>	12 AH Lithium-Ion
<b>Colors</b>	Red, Blue
<b>Freewheel Mode</b>	Yes
<b>Electro-Mechanical Brakes</b>	Yes
<b>Charger</b>	2A, Off Board

<b>Controller</b>	40 Amp
<b>Length</b>	38"
<b>Width</b>	24"
<b>Height (ground to top of back of seat)</b>	34.25"
<b>Ground to Top of Seat</b>	21"
<b>Weight of Unit Assembled</b>	74 lbs.
<b>Weight of Base</b>	46 lbs.
<b>Weight of Seat</b>	22 lbs.
<b>Weight of Battery</b>	6 lbs.
<b>Drive Wheels</b>	7.75"
<b>Standard Seat Width x Depth</b>	18" x 17"
<b>Back Height - No Headrest</b>	15.5"
<b>Color</b>	Black
<b>Warranty Frame</b>	5 Years
<b>Warranty Drive Train</b>	13 Months
<b>Warranty Electronics</b>	13 Months
<b>Warranty Batteries</b>	12 Months