

**COMPASS™ HD****\$6,599.00**

The Golden Compass HD offers you a stable, center-wheel-drive platform and an articulating chassis for stability over tough terrain. A bigger, beefier version of the popular Compass Sport, the Compass HD has a 450 lb. weight capacity, yet a modest 20-inch turning radius.

SKU: GP620**Categories:** [Center Wheel Power Chairs](#),
[Power Wheelchairs](#)**ADDITIONAL INFORMATION**

Length	40"
Width	25"
Front Tire	N/A
Rear tire	N/A
Center Tire	10"
Casters	6"
Weight Capacity	450 lbs
Maximum Speed	4 mph
Ground Clearance	3.5"
Operating Range	12 miles
Turning Radius	20"



Batteries	NF22
Weight w batteries	226 lbs
Front Half Weight	N/A
Rear Half Weight	N/A
Seat weight w arms	51 lbs
Heaviest Piece	101 lbs
Base Weight	101 lbs
Features	Adjustable Armrests, Adjustable Headrest, Adjustable Seat Posts, Center Wheel Drive, Electronic Speed Control, Freewheel Mode, Programmable Joystick, Slide Seat Bracket
Warranty on Electrical Parts	13 months
Warranty Batteries	1 year