

**COMPASS™ HD****\$6,599.00**

The Golden Compass HD offers you a stable, center-wheel-drive platform and an articulating chassis for stability over tough terrain. A bigger, beefier version of the popular Compass Sport, the Compass HD has a 450 lb. weight capacity, yet a modest 20-inch turning radius.

---

**SKU:** GP620**Categories:** [Center Wheel Power Chairs](#),  
[Power Wheelchairs](#)**ADDITIONAL INFORMATION**

|                         |                         |
|-------------------------|-------------------------|
| <b>Length</b>           | 40"                     |
| <b>Width</b>            | 25"                     |
| <b>Front Tire</b>       | N/A                     |
| <b>Rear tire</b>        | N/A                     |
| <b>Center Tire</b>      | 10"                     |
| <b>Casters</b>          | 6"                      |
| <b>Weight Capacity</b>  | <a href="#">450 lbs</a> |
| <b>Maximum Speed</b>    | 4 mph                   |
| <b>Ground Clearance</b> | 3.5"                    |
| <b>Operating Range</b>  | 12 miles                |
| <b>Turning Radius</b>   | 20"                     |



|                                     |  |
|-------------------------------------|--|
| <b>Batteries</b>                    | NF22   |
| <b>Weight w batteries</b>           | 226 lbs  |
| <b>Front Half Weight</b>            | N/A  |
| <b>Rear Half Weight</b>             | N/A  |
| <b>Seat weight w arms</b>           | 51 lbs   |
| <b>Heaviest Piece</b>               | 101 lbs  |
| <b>Base Weight</b>                  | 101 lbs  |
| <b>Features</b>                     | Adjustable Armrests, Adjustable Headrest, Adjustable Seat Posts, Center Wheel Drive, Electronic Speed Control, Freewheel Mode, Programmable Joystick, Slide Seat Bracket |
| <b>Warranty on Electrical Parts</b> | 13 months  |
| <b>Warranty Batteries</b>           | 1 year   |