



COMPASS™ HD

\$6,599.00

The Golden Compass HD offers you a stable, center-wheel-drive platform and an articulating chassis for stability over tough terrain. A bigger, beefier version of the popular Compass Sport, the Compass HD has a 450 lb. weight capacity, yet a modest 20-inch turning radius.

FDA Class II Medical Devices designed by Golden are intended to assist persons with chronic disablements.

SKU: GP620

Categories: [Center Wheel Power Chairs](#),
[Power Wheelchairs](#)

ADDITIONAL INFORMATION

Medicare Code	K0825
User Weight Capacity	450 lbs
Drive Wheels	Center
Maximum Speed	4 mph
Operating Range	12 miles
Ground Clearance	3.5"
Turning Radius	20"
Type Batteries	NF22
Standard Colors	Red, Blue



Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	75 Amp
Length	40"
Width	25"
Height (ground to top of back of seat)	38.25" - 42.75"
Height (ground to top of back of headrest)	47"-53"
Ground to Top of Seat	21.25" - 24.25"
Number Seat Height Positions	5
Size of Increments	0.5"
Weight of Unit (assembled)	226 lbs
Base Weight	101 lbs
Seat with Arms Weight	51 lbs
Batteries Weight	37 lbs ea
Tires (standard)	Flat Free / Non Marking Solid
Center Tire	10"
Casters	6"
Standard Seat	Captain's Seat
Back Height (no headrest)	19"
Back Height with Headrest	26.5"-30"
Seat Width x Depth	18"-24"x18"-20"
Seat Color	Black Vinyl
Optional Seat	Captain's Chair
Optional Seat 2 Back Height (no headrest)	19"



Optional Seat 2 Back Height (with headrest)	26.5"-30"
Optional Seat Width x Depth	18"-24"x18"-20"
Optional Seat Color	Black Vinyl
Footrest Height Adjustable	Yes
Footrest Angle Adjustable	Yes
Warranty Frame	5 Years
Warranty Drive Train	13 months
Warranty Electronics	13 months
Warranty Batteries	1 year