

FDA Class II Medical Device



COMPASS™ SPORT

\$4,999.00

The Golden Compass™ Sport combines sleek style and superior maneuverability with a smooth, intuitive ride, thanks to the Dynamic LiNX Smart Controller and patented springless articulating chassis—plus customizable comfort with multiple seat widths and a luxury Captain's seat.

SKU: GP605

Categories: [Center Wheel Power Chairs](#),
[Power Wheelchairs](#)

ADDITIONAL INFORMATION

Medicare Code	K0823
User Weight Capacity	300 lbs
Drive Wheels	Center
Maximum Speed	4 mph
Operating Range	24 miles
Ground Clearance	3.25"
Turning Radius	24"
Type Batteries	U1
Optional Batteries	50 AH
Standard Colors	Red, Blue
Freewheel Mode	Yes



Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	40 Amp
Length	39.5"
Width	24"
Height (ground to top of back of seat)	40" - 42"
Height (ground to top of back of headrest)	47"-53"
Ground to Top of Seat	19" - 21"
Weight of Unit (assembled)	179 lbs
Base Weight	97 lbs
Seat with Arms Weight	40 lbs
Center Tire	10"
Casters	6"
Standard Seat	Captain's Seat
Back Height (no headrest)	23"
Back Height with Headrest	26.5"-30"
Seat Width x Depth	16"-22"x16"-20"
Seat Color	Black Vinyl
Optional Seat	Captain's Chair
Optional Seat 2 Back Height (no headrest)	19"
Optional Seat 2 Back Height (with headrest)	26.5"-30"
Optional Seat Width x Depth	16"-22"x16"-20"
Optional Seat Color	Black Vinyl
Footrest Height Adjustable	Yes
Footrest Angle Adjustable	Yes



Warranty Frame	Lifetime
Warranty Drive Train	13 months
Warranty Electronics	13 months
Warranty Batteries	1 year