

**COMPASS™ SPORT****\$4,999.00**

The Golden Compass™ Sport features a stylish, sleek design built for maneuverability both indoors and out with durable, solid color injected plastic shrouds in red or blue for a sporty new look. It offers a range of comfortable seat sizes available in 16", 18" and 20" standard widths and an upgraded 22" wide seat.

FDA Class II Medical Devices designed by Golden are intended to assist persons with chronic disablements.

---

**SKU:** GP605**Categories:** [Center Wheel Power Chairs](#),  
[Power Wheelchairs](#)**ADDITIONAL INFORMATION**

<b>Medicare Code</b>	K0823
<b>User Weight Capacity</b>	<a href="#">300 lbs</a>
<b>Drive Wheels</b>	Center
<b>Maximum Speed</b>	4 mph
<b>Operating Range</b>	24 miles
<b>Ground Clearance</b>	3.25"
<b>Turning Radius</b>	24"
<b>Type Batteries</b>	U1



<b>Optional Batteries</b>	50 AH
<b>Standard Colors</b>	Red, Blue
<b>Freewheel Mode</b>	Yes
<b>Electronic Speed Control</b>	Yes
<b>Electro-Mechanical Brakes</b>	Yes
<b>Charger</b>	Off Board
<b>Controller</b>	40 Amp
<b>Length</b>	39.5"
<b>Width</b>	24"
<b>Height (ground to top of back of seat)</b>	40" - 42"
<b>Height (ground to top of back of headrest)</b>	47"-53"
<b>Ground to Top of Seat</b>	19" - 21"
<b>Weight of Unit (assembled)</b>	179 lbs
<b>Base Weight</b>	97 lbs
<b>Seat with Arms Weight</b>	40 lbs
<b>Center Tire</b>	10"
<b>Casters</b>	6"
<b>Standard Seat</b>	Captain's Seat
<b>Back Height (no headrest)</b>	23"
<b>Back Height with Headrest</b>	26.5"-30"
<b>Seat Width x Depth</b>	16"-22"x16"-20"
<b>Seat Color</b>	Black Vinyl
<b>Optional Seat</b>	Captain's Chair
<b>Optional Seat 2 Back Height (no headrest)</b>	19"
<b>Optional Seat 2 Back Height (with headrest)</b>	26.5"-30"
<b>Optional Seat Width x Depth</b>	16"-22"x16"-20"



<b>Optional Seat Color</b>	Black Vinyl
<b>Footrest Height Adjustable</b>	Yes
<b>Footrest Angle Adjustable</b>	Yes
<b>Warranty Frame</b>	Lifetime
<b>Warranty Drive Train</b>	13 months
<b>Warranty Electronics</b>	13 months
<b>Warranty Batteries</b>	1 year