

FDA Class II Medical Device



## COMPASS™ SPORT

**\$4,999.00**

The Golden Compass™ Sport combines sleek style and superior maneuverability with a smooth, intuitive ride, thanks to the Dynamic LiNX Smart Controller and patented springless articulating chassis—plus customizable comfort with multiple seat widths and a luxury Captain's seat.

---

**SKU:** GP605

**Categories:** [Center Wheel Power Chairs](#),  
[Power Wheelchairs](#)

## ADDITIONAL INFORMATION

<b>Medicare Code</b>	K0823
<b>User Weight Capacity</b>	<a href="#">300 lbs</a>
<b>Drive Wheels</b>	Center
<b>Maximum Speed</b>	4 mph
<b>Operating Range</b>	24 miles
<b>Ground Clearance</b>	3.25"
<b>Turning Radius</b>	24"
<b>Type Batteries</b>	2-U1/ up to 50AH
<b>Optional Batteries</b>	50 AH
<b>Standard Colors</b>	Red, Blue
<b>Freewheel Mode</b>	Yes



<b>Electronic Speed Control</b>	Yes
<b>Electro-Mechanical Brakes</b>	Yes
<b>Charger</b>	Off Board
<b>Controller</b>	40 Amp
<b>Length</b>	39.5"
<b>Width</b>	24"
<b>Height (ground to top of back of seat)</b>	40" - 42"
<b>Height (ground to top of back of headrest)</b>	47"-53"
<b>Ground to Top of Seat</b>	19" - 21"
<b>Weight of Unit (assembled)</b>	179 lbs
<b>Base Weight</b>	97 lbs
<b>Seat with Arms Weight</b>	40 lbs
<b>Casters</b>	6"
<b>Standard Seat</b>	Captain's Seat
<b>Back Height (no headrest)</b>	23"
<b>Back Height with Headrest</b>	30.5" - 34.5"
<b>Seat Width x Depth</b>	18" x 18"/20"
<b>Seat Color</b>	Black Vinyl
<b>Optional Seat</b>	Captain's Chair
<b>Optional Seat 2 Back Height (no headrest)</b>	19"
<b>Optional Seat 2 Back Height (with headrest)</b>	26.5"-30"
<b>Optional Seat Width x Depth</b>	16", 18", 20", or 22" x 16", 18" or 20"
<b>Optional Seat Color</b>	Black Vinyl
<b>Footrest Height Adjustable</b>	Yes
<b>Footrest Angle Adjustable</b>	Yes
<b>Warranty Frame</b>	Lifetime

<b>Warranty Drive Train</b>	13 months
<b>Warranty Electronics</b>	13 months
<b>Warranty Batteries</b>	1 year