



COMPASS™ SPORT

\$4,999.00

The Golden Compass™ Sport features a stylish, sleek design built for maneuverability both indoors and out with durable, solid color injected plastic shrouds in red or blue for a sporty new look. It offers a range of comfortable seat sizes available in 16", 18" and 20" standard widths and an upgraded 22" wide seat.

SKU: GP605

Categories: [Center Wheel Power Chairs](#),
[Power Wheelchairs](#)

ADDITIONAL INFORMATION

Length	39.5"
Width	24"
Front Tire	N/A
Rear tire	N/A
Center Tire	10"
Casters	6"
Weight Capacity	300 lbs
Maximum Speed	4 mph
Ground Clearance	3.25"
Operating Range	24 miles



Turning Radius	21"
Batteries	U1
Weight w batteries	179 lbs
Front Half Weight	N/A
Rear Half Weight	N/A
Seat weight w arms	40 lbs
Heaviest Piece	97 lbs
Base Weight	97 lbs
Warranty Batteries	1 year
Warranty on Electrical Parts	13 months
Warranty Drive Train	13 months
Features	Adjustable Armrests, Adjustable Headrest, Center Wheel Drive, Electronic Speed Control, Footrest Extension, Freewheel Mode, Programmable Joystick, Slide Seat Bracket