

FDA Class II Medical Device



## CRICKET LT CARBON

GO LIGHT. FOLD FAST. SPEND LESS. At just 25.8 lbs., the Golden Cricket LT Carbon folds sideways for easy transport while supporting up to 300 lbs.

---

**SKU:** GP309

**Categories:** [Cricket Power Chairs](#), [Folding, Power Wheelchairs](#), [Travel Power Chairs](#)

**Tags:** [New](#)

## ADDITIONAL INFORMATION

<b>Medicare Code</b>	K0899
<b>User Weight Capacity</b>	<a href="#">300 lbs</a>
<b>Drive Wheels</b>	Rear
<b>Maximum Speed</b>	3.7 mph
<b>Operating Range</b>	9.5 miles
<b>Ground Clearance</b>	1.25"
<b>Turning Radius</b>	29.5"
<b>Type Batteries</b>	10AH 24V Lithium-Ion
<b>Color</b>	Black Carbon Fiber
<b>Freewheel Mode</b>	Yes
<b>Electronic Speed Control</b>	Yes
<b>Electro-Mechanical Brakes</b>	Yes
<b>Charger</b>	Off Board

<b>Controller</b>	35 Amp
<b>Length</b>	36"
<b>Width</b>	24"
<b>Height (ground to top of back of seat)</b>	33"
<b>Ground to Top of Seat</b>	19"
<b>Number of Seat Height Positions</b>	N/A
<b>Size of Increments</b>	N/A
<b>Weight of Unit (assembled)</b>	31.3 lbs
<b>Battery Weight</b>	2 lbs
<b>Tires</b>	Flat Free
<b>Drive Wheels (Rear)</b>	11"
<b>Casters (Front)</b>	6.5"
<b>Rear Anti-Tip Wheels</b>	2"
<b>Standard Seat Back Height - no headrest</b>	16.4"
<b>Seat Width x Depth</b>	18" x 18"
<b>Seat Color</b>	Black
<b>Footrest Height Adjustable</b>	N/A
<b>Footrest Angle Adjustable</b>	N/A
<b>Warranty Frame</b>	5 Years
<b>Warranty Drive Train</b>	13 months
<b>Warranty Electronics</b>	13 months
<b>Warranty Battery</b>	1 Year