



CRICKET

\$3,799.00

The Cricket Power Wheelchair is stylish, agile, and extremely portable. Crafted with an incredibly strong yet lightweight 100% carbon fiber frame, the Cricket weighs as light as 33 lbs. making travel a breeze! The slim sporty design is perfect for navigating indoors and out while sightseeing, shopping, or dining.

SKU: GP302

Categories: [Foldable Power Wheelchairs](#), [Power Wheelchairs](#), [Travel Power Chairs](#)

ADDITIONAL INFORMATION

Medicare Code	Not Coded
User Weight Capacity	300 lbs
Drive Wheels	Rear
Maximum Speed	3.7 mph
Operating Range	9.3 miles
Ground Clearance	2.6"
Turning Radius	22"
Type Batteries	12AH 24V Lithium-Ion
Optional Batteries	N/A
Standard Colors	Black



Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	35 Amp
Length	36.6"
Width	22.5"
Height (ground to top of back of seat)	36.5"
Height (ground to top of back of headrest)	N/A
Ground to Top of Seat	19.3"
Ground to Top of Deck	N/A
Deck to Top of Seat	N/A
Weight of Unit (assembled)	39.3 lbs
Tires (standard)	Flat Free / Non Marking Solid
Rear Tires	8.5"
Casters	7"
Rear Anti-Tippers	1.6"
Standard Seat	Mesh Back/ Removable Seat Cushion
Back Height (no headrest)	19.3"
Seat Width x Depth	17.3"x16.5"
Seat Color	Black
Warranty Frame	5 Years
Warranty Drive Train	13 months
Warranty Electronics	13 months
Warranty Batteries	1 year