

FDA Class II Medical Device



## CRICKET

**\$3,799.00**

The Cricket Power Wheelchair is stylish, agile, and extremely portable. Crafted with an incredibly strong yet lightweight 100% carbon fiber frame, the Cricket weighs as light as 33 lbs. making travel a breeze! The slim sporty design is perfect for navigating indoors and out while sightseeing, shopping, or dining.

---

**SKU:** GP302

**Categories:** [Favorites](#), [Foldable Power Wheelchairs](#), [Power Wheelchairs](#), [Travel Power Chairs](#)

**Tags:** [Best Seller](#)

## ADDITIONAL INFORMATION

<b>Medicare Code</b>	Not Coded
<b>User Weight Capacity</b>	<a href="#">300 lbs</a>
<b>Drive Wheels</b>	Rear
<b>Maximum Speed</b>	3.7 mph
<b>Operating Range</b>	9.3 miles
<b>Ground Clearance</b>	2.6"
<b>Turning Radius</b>	22"
<b>Type Batteries</b>	12AH 24V Lithium-Ion
<b>Optional Batteries</b>	N/A
<b>Standard Colors</b>	Black



<b>Freewheel Mode</b>	Yes
<b>Electronic Speed Control</b>	Yes
<b>Electro-Mechanical Brakes</b>	Yes
<b>Charger</b>	Off Board
<b>Controller</b>	35 Amp
<b>Length</b>	36.6"
<b>Width</b>	22.5"
<b>Height (ground to top of back of seat)</b>	36.5"
<b>Height (ground to top of back of headrest)</b>	N/A
<b>Ground to Top of Seat</b>	19.3"
<b>Ground to Top of Deck</b>	N/A
<b>Deck to Top of Seat</b>	N/A
<b>Weight of Unit (assembled)</b>	39.3 lbs
<b>Tires (standard)</b>	Flat Free / Non Marking Solid
<b>Rear Tires</b>	8.5"
<b>Casters</b>	7"
<b>Rear Anti-Tippers</b>	1.6"
<b>Standard Seat</b>	Mesh Back/ Removable Seat Cushion
<b>Back Height (no headrest)</b>	19.3"
<b>Seat Width x Depth</b>	17.3"x16.5"
<b>Seat Color</b>	Black
<b>Warranty Frame</b>	5 Years
<b>Warranty Drive Train</b>	13 months
<b>Warranty Electronics</b>	13 months
<b>Warranty Batteries</b>	1 year