



LITERIDER ENVY LT

\$2,299.00

The LiteRider Envy LT Power Wheelchair is both stylish and durable with a tubular design that is unique to the Golden LiteRider family! The Envy LT features a compact PG Drives Technology nVR2 controller and well-styled joystick designed to make everyday tasks easier such as driving up to a table or maneuvering through a doorway.

SKU: GP161

Categories: [Compact Power Chairs](#), [Power Wheelchairs](#), [Travel Power Chairs](#)

ADDITIONAL INFORMATION

Length	33.25"
Width	22.5"
Front Tire	N/A
Rear tire	8"
Center Tire	N/A
Casters	6"
Weight Capacity	250
Maximum Speed	3.5 MPH
Ground Clearance	2.5"
Operating range	7.25 MILES



Turning radius	28.5"
Batteries	18 AH
Weight w batteries	113 LBS
Front Half Weight	34.5 LBS
Rear Half Weight	30.5 LBS
Seat weight w arms	23 LBS
Heaviest Piece	34.5 LBS
Base Weight	N/A
Warranty Drive Train	13 Months
Warranty Batteries	1 Year
Features	Adjustable Seat Posts, Electronic Speed Control, Footrest Extension, Freewheel Mode