

FDA Class II Medical Device



LITERIDER ENVY

\$2,239.00

Designed with the end-user and their caregivers in mind, the LiteRider Envoy can be disassembled in just seconds to easily fit in the trunk of virtually any car. Available in 6 colors!

FDA Class II Medical Devices designed by Golden are intended to assist persons with chronic disabilities.

SKU: GP162

Categories: [Compact Power Chairs](#), [Power Wheelchairs](#), [Travel Power Chairs](#)

Tags: [Best Seller](#)

ADDITIONAL INFORMATION

Medicare Code	Not Coded
User Weight Capacity	300 lbs
Drive Wheels	Rear
Maximum Speed	3.5 mph
Operating Range	15.5 miles
Ground Clearance	2.5"
Turning Radius	28.5"
Type Batteries	2-22 AH
Standard Colors	Red, Blue, Silver, Green, Yellow, Orange
Freewheel Mode	Yes



Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	40 Amp
Length	36"
Width	22.5"
Height (ground to top of back of seat)	34" - 35"
Ground to Top of Seat	21.5" - 22.5"
Weight of Unit (assembled)	123 lbs
Front Half Weight	37 lbs
Rear Half Weight	36 lbs
Seat with Arms Weight	23 lbs
Batteries Weight	27 lbs
Tires (standard)	Flat Free / Non Marking Solid
Rear Tires	9"
Casters	6"
Rear Anti-Tippers	3"
Standard Seat	Highback
Back Height (no headrest)	16"
Seat Width x Depth	18" x 16"
Seat Color	Black Vinyl
Optional Seat	Stadium Style
Optional Seat 2 Back Height (no headrest)	16"
Optional Seat 2 Width x Depth	20" x 18"
Optional Seat Color	Black Vinyl
Footrest Height Adjustable	Yes
Footrest Angle Adjustable	Yes



Warranty Frame	5 Years
Warranty Drive Train	13 months
Warranty Electronics	13 months
Warranty Batteries	1 year