

FDA Class II Medical Device



LITERIDER ENVY

\$2,239.00

Designed with the end-user and their caregivers in mind, the LiteRider Envy can be disassembled in just seconds to easily fit in the trunk of virtually any car. Available in 6 colors!

SKU: GP162

Categories: [Compact Power Chairs](#), [Power Wheelchairs](#), [Travel Power Chairs](#)

Tags: [Best Seller](#)

ADDITIONAL INFORMATION

Length	36"
Width	22.5"
Front Tire	N/A
Rear tire	9"
Center Tire	N/A
Casters	6"
Weight Capacity	300 lbs
Maximum Speed	3.5 mph
Ground Clearance	2.5"
Operating Range	15.5 miles
Turning Radius	28.5"
Batteries	22 AH
Weight w batteries	123 lbs



Front Half Weight	37 lbs
Rear Half Weight	36 lbs
Seat weight w arms	23 lbs
Heaviest Piece	37 lbs
Base Weight	N/A
Warranty Batteries	1 year
Warranty Drive Train	13 months
Warranty on Electrical Parts	13 months
Features	Adjustable Armrests, Adjustable Headrest, Adjustable Seat Posts, Center Wheel Drive, Electronic Speed Control, Footrest Extension, Freewheel Mode, Premium Seat Option