

FDA Class II Medical Device



STRIDE

\$2,499.00

The GP301 Stride Power Wheelchair will take you where you want to go in comfort and style! The Stride is stylish with a two-tone grey and black seat with removable cushion and a tall back for long-lasting comfort while on the go.

SKU: GP301

Categories: [Compact Power Chairs](#), [Foldable Power Wheelchairs](#), [Power Wheelchairs](#), [Travel Power Chairs](#)

ADDITIONAL INFORMATION

Height (ground to top of back of seat)	39"
User Weight Capacity	300 lbs
Height (ground to top of back of headrest)	N/A
Maximum Speed	3.7 mph
Operating Range	9.3 miles
Ground Clearance	1.6"
Turning Radius	47.2"
Type Batteries	10AH 24V Lithium-Ion
Standard Color	Black
Freewheel Mode	Yes



Electronic Speed Control	Yes
Charger	Off Board
Controller	50 Amp
Length	42.5"
Ground to Top of Seat	22"
Ground to Top of Deck	N/A
Deck to Top of Seat	N/A
Number Seat Height Positions	N/A
Size of Increments	N/A
Front Axle to Rear Axle	N/A
Front Half Weight	N/A
Rear Half Weight	N/A
Base Weight	N/A
Seat with Arms Weight	N/A
Batteries Weight	3.8 lbs
Tires (standard)	Flat Free / Non Marking Solid
Rear Tires	12"
Casters	8"
Standard Seat	Sling Back/ Removable Seat Cushion
Back Height (no headrest)	16"
Back Height with Headrest	N/A
Seat Width x Depth	17"x18"
Seat Color	Black
Warranty Frame	5 Years
Warranty Drive Train	13 months
Warranty Electronics	13 months
Warranty Batteries	1 year