



STRIDE

\$2,799.00

The GP301 Stride Power Wheelchair will take you where you want to go in comfort and style! The Stride is stylish with a two-tone grey and black seat with removable cushion and a tall back for long-lasting comfort while on the go.

SKU: GP301

Categories: [Compact Power Chairs](#), [Foldable Power Wheelchairs](#), [Power Wheelchairs](#), [Travel Power Chairs](#)

Tags: [New](#)

ADDITIONAL INFORMATION

Length	42.5"
Width	24.4"
Front Tire	8"
Rear tire	12"
Weight Capacity	300 lbs
Maximum Speed	3.7 mph
Ground Clearance	1.6"
Operating Range	9.3 miles
Turning Radius	47.2"
Batteries	10AH 24V Lithium-Ion
Weight w batteries	58 lbs

Warranty on Electrical Parts

13 months

Features

50 Amp Controller, Airline Approved, Aluminum Frame, Comfort Spring Suspension, Easy Storage, Joystick with USB Charger, Removable Padded Seat, Under-Seat & Storage & Back Seat Pocket

Warranty Batteries

1 year

Rear Anti-Tip Wheels

2"

Warranty Drive Train

13 Months

Standard Color

Black