



Standard Colors:



Crimson Red



Azure Blue

NEW FOR 2024!

- NEW High Back Stadium Seat with Quilted Diamond Stitch
- NEW Azure Blue & Crimson Red Color Shrouds

Companion 3-Wheel (Mid-Size) GC240

GOLDEN
Companion

The Companion scooter provides great style with high performance and user-friendly disassembly. Features include front and rear lighting package with side reflectors; LED battery gauge display; and wraparound Delta tiller with lock mechanism. The Companion offers great value in a full-time portable scooter.

Standard Features:

- Under Mat Troubleshooting Guide
- Delta Tiller with Auto Lock Mechanism
- Easy Reach Tiller Adjustment
- Adjustable Sliding Seat
- Adjustable LED Headlight
- U-1 Batteries



GOLDEN Companion

GC240


Companion 3-Wheel Specifications

Model	GC240
Medicare Code	K0807
Weight Capacity	350 lbs
Drive Wheels	Rear
Maximum Speed	5 mph
Operating Range ¹	15 miles
Ground Clearance Under Deck	3.75"
Ground Clearance to Transaxle/Motor	2.5"
Turning Radius	37"
Type Batteries / Air Travel Safe	2-U1 / Yes
Colors	Red, Blue
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Energy Efficient - Auto Turn Off (Sleep Mode)	Yes
Charger	Portable 3.5 Amp Off Board
Controller	90 Amp Dynamic R Series
Dimensions:	
Length	42.25"
Width	24"
Ground to Top of Seat	22.25" - 24.30"
Ground to Top of Deck	6.25"
Deck to Top of Seat	16.25" - 18.25"
Number of Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	30"
Weight of Unit: (assembled)	170.8 lbs
Tires:	Flat Free / Non Marking Solid
Front	9"
Rear	10.75"
Front & Rear Anti-Tippers	3"
Standard Seat:	Quilted High Back
Slide Seat	Standard
Back Height - no headrest	18"
Width x Depth	18" x 16"
Color	Black
Transportation Seat Lock	Yes, Self Locking
Seat Rotation	90 Degree either direction 360 Degrees
Optional Seat:	Stadium Style
Back Height - no headrest	16"
Width x Depth	20" x 18"
Color	Black
Seat Options:	
Pneumatic Rise Seat	Yes
Warranty:	
Frame	Lifetime
Drive Train	2 Years
Electronics ²	2 Years
Batteries	1 Year

MID-SIZE 3-WHEEL SCOOTER




User-friendly auto-tiller adjustment



Infinite adjustable tiller with ergonomic design



Front and rear lights and side reflectors for safety



Trouble Shooting Guide conveniently located under mat



Disassembles easily!

¹ Battery range will vary due to rider weight, drive surface and terrain.

² Electronics warranty excludes batteries.



Phone: 800-624-6374

Fax: 800-628-5165

www.goldentech.com

Find us on Facebook

www.facebook.com/goldentechnologies