




Standard Color



Black

AS
LIGHT AS
52.6
lbs.!*

Stride GP301



The all-new Stride Power Wheelchair will take you where you want to go in comfort and style! The Stride is stylish with a two-tone grey and black seat with removable cushion and a tall back for long-lasting comfort while on the go. Modern-looking front and back comfort spring suspension add to the smooth ride. The sporty design is perfect for navigating indoors and out while sightseeing, shopping, or dining. This portable folding power chair is lightweight with an aluminum frame and stows easily in the trunk of your car. Weighing as light as 52.6 lbs.*, the Stride is the perfect choice to take on road trips, cruise ships and airplanes. The Stride features a well-styled joystick designed to make everyday tasks easier such as driving up to a table or maneuvering through a doorway. It has an impressive 23.6" turning radius and offers a 300 lb. weight capacity, maximum speed of 3.7 mph, and an operating range of up to 9.3 miles on airline-approved 10AH 24V lithium-ion batteries. It also offers convenient offboard charging and folds effortlessly for transport.

Standard Features:

- Aluminum Frame
- 300 lb. Weight Capacity
- Removable Compact Joystick with USB Charger
- 50 Amp Controller
- Comfort Spring Suspension
- Removable Padded Seat
- Convenient Under-Seat Storage & Back Seat Pocket
- Offboard Charging Capability
- Airline-approved Lithium-ion Battery
- Easy One-Step Folding
- Stands on End for Storage





FOLDABLE POWER CHAIR

Model Number	GP301
Medicare Code	Not Coded
Weight Capacity	300 lbs.
Drive Wheels	Rear
Maximum Speed	3.7 mph
Operating Range ¹	9.3 miles
Ground Clearance	1.6"
Turning Radius	23.6"
Type Batteries	10AH 24V Lithium-Ion
Color	Black
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	2 Amp, Off Board
Controller	50 Amp
Dimensions:	
Length	42.5"
Width	24.4"
Height (ground to top of back of seat)	38.6"
Ground to Top of Seat	22"
Number of Seat Height Positions	N/A
Size of Increments	N/A
Weight of Unit: (assembled)	58.0 lbs.
Weight of Unit: (minus parts)*	52.6 lbs.
Front Half	N/A
Rear Half	N/A
Seat w/Arms	N/A
Batteries	3.8 lbs
Tires:	
Drive Wheels (Rear)	12"
Casters (Front)	8"
Rear Anti-Tip Wheels	2"
Standard Seat:	
	Sling Back with Removable Seat Cushion
Back Height - no headrest	16"
Width x Depth	17" x 18"
Color	Black
Footrest:	
Height adjustable	N/A
Angle adjustable	N/A
Warranty:	
Frame	5 Years
Drive Train	13 Months
Electronics ²	13 Months
Batteries	12 Months

¹ Battery range tested at 200 lbs. and will vary due to rider weight, drive surface, terrain and battery type.

² Electronics warranty excludes batteries.

* Excludes battery, joystick, cushion and bag.

FDA Class II Medical Device

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

AS LIGHT AS 52.6 lbs.!*



Portable & high-performing with an aluminum frame!



FOLD & STOW!
Compact Dimensions
42.5" x 24.4" x 15.7"

East One-Step Folding with Velcro Closure Strap for Transport

Stands On End for Storage



Comfort Spring Suspension for a Smooth Ride



Removable Padded Seat



Phone: 800-624-6374

Fax: 800-628-5165

www.goldentech.com



www.facebook.com/goldentechnologies

SHT-GP301-031524