

Buzzaround 3 wheel GB121-SHZ 4 wheel GB124-SHZ



The next generation of the popular Buzzaround XLSHD features compact portability with a heavy-duty weight capacity! The XLSHD features a robust 325 lb. weight capacity and our Comfort Spring Suspension for a smoother ride. It offers improved standard features like a larger wraparound Delta tiller, higher ground clearance of 4", and an operating range of up to 9 miles on the standard 12 AH batteries. Upgrade with the new XR Extended Range Battery option for even more mileage! The Buzzaround XLSHD also offers a new stylish, contoured stadium style seat for more comfort and increased seating tolerance. This scooter disassembles in seconds with just four easy steps making transporting it from one destination to another a breeze!



Standard Colors:





(Included in Box)

NEW Sporty Optional Colors:



Orchid Pink

(Matte)

Desert Tan

Features:

- 325 lb. Weight Capacity
- Comfort Spring Suspension
- Manually Adjustable Preload on Rear Suspension
- Improved Drivability
- Higher Ground Clearance
- LED Angle-adjustable Headlight
- USB Port on the Tiller
- Phone & Cup Holders
- Quick & Easy Disassembly

Upgrade with NEW **XR Extended Range Battery Option!**

(Packaged separately as an accessory but ships with the unit.)



LIGHTWEIGHT HEAVY-DUTY PORTABLE SCOOTER

Item	GB121-SHZ	GB124-SHZ
Description	XLSHD	XLSHD
Medicare Code	K0801	K0801
Weight Capacity	325 lbs.	325 lbs.
Drive Wheels	Rear	Rear
Maximum Speed	3.95 mph	3.7 mph
Operating Range up to ¹	6.3 miles	6.98 miles
Under Deck Clearance	4"	4.5"
Turning Radius	36.5"	49.5"
Type Batteries	2-12 AH	2-12 AH
Colors	Red, Blue	Red, Blue
Freewheel Mode	Yes	Yes
Electronic Speed Control	Yes	Yes
Electro-Mechanical Brakes	Yes	Yes
Charger	Off-board	Off-board
Controller	50 Amp	50 Amp
Dimensions		
Length	41"	40.5"
Width	22.25"	21.5"
Height (ground to top of back of seat)	36.25"-37.75"	33.5"-35.5"
Ground to Top of Seat	21"-22.5"	18.5"-20.5"
Ground to Top of Deck	5.5"	4.5"
Deck to Top of Seat	15.5"-17"	14"-16"
Number of Seat Height Posi	itions 3	3
Size of Increments	.75"	1"
Front Axle to Rear Axle	31.5"	30.5"
Weight of Unit: (assembled) ³	106.5 lbs.	103.5 lbs.
Front Half	35 lbs.	40 lbs.
Rear Half	28 lbs.	28 lbs.
Seat w/Arms	21 lbs.	19.5 lbs.
Batteries (pack with 2- bat	teries) 16 lbs	16 lbs. lbs.
Tires:	Solid, low profile	Solid, low profile
Front	200 x 50 mm	190 x 54 mm
Rear	198 x 70 mm	198 x 70 mm
Rear Anti-Tip Wheels	2"	2"
Standard Seat	Stadium Style	Stadium Style
Seat Back Height - no headrest	15.25"	15"
Width x Depth	16.5" x 15.13"	16.5" x 15"
Color	Black Vinyl	Black Vinyl
Warranty:	,	
Frame	Lifetime Warranty	Lifetime Warranty
	,	
Drive Train	2 Year Warranty	2 Year Warranty
Drive Train Electronics ²	2 Year Warranty 2 Year Warranty	2 Year Warranty 2 Year Warranty

^{1.} Battery range at 200 lbs. and will vary due to rider weight, drive surface, terrain, and battery type.

^{2.} Electronics warranty excludes batteries.





Disassembles into 5 Pieces!



Comfort Spring Suspension Manually adjustable preload on rear suspension!



Phone: 800-624-6374 Fax: 800-628-5165 www.goldentech.com



Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

SHT-GB121-GB124-SHZ-101624