





Standard Fabrics:



Calypso (TCP)




Hazelnut (THA)



Sterling (TST)

The Tucker pattern has a modern, updated look with rich color options. It features our exclusive SofTouch™ with Stain Defense coating.



Umber (SUM)*

The Sinda polyurethane pattern is soft and supple. It was designed to be breathable yet highly durable.

5-Year Sinda Warranty!

*No upcharge for Sinda on PR401.

Crossover fabrics are unavailable.

Cambridge PR401



Comfort Zone[™]
by GOLDEN

Enjoy the simplicity of comfort in the Cambridge by Golden. This traditional style power lift recliner is available in three sizes offering a unique fit for hours of comfort. The plush biscuit back and luxurious chaise pad provide supportive comfort for your body while the 1-Zone Comfort Zone Positioning System glides your body to an ergonomic recline position.

Standard Features:

- Classic Design
- Blissful Comfort with Ultra Soft Back
- Simple Two-Button Hand Control
- SmartTek® Diagnostic System



DeLUNA Series PR401



Model Number	PR401-SME	PR401-MLA
Size	Small/Med	Med/Large
Class II Medical Device	YES	YES
Weight Capacity	375 Lbs.	375 Lbs.
Positioning	3	3
Comfort Zone(s)	1	1
Twilight Technology	NO	NO
ZG+	NO	NO
Chaise Seat	YES	YES
Pocketed Coil Seat	NO	NO
Back Type	Biscuit	Biscuit
Removable Back Type	Quick Release	Quick Release
Head Rest Cover	Standard	Standard
Arm Rest Covers	Standard	Standard
Overall Width	33"	35"
Overall Height	41.5"	44.5"
Seat Width	20"	22"
Floor to Top of Seat	19"	22"
Seat Depth	19"	22"
Seat to Top of Back	27"	31"
Seat to Top of Arm	6.5"	7"
Distance Required from Wall Reclined	18"	19"
Extension reclined from Back of Seat to End of Foot Rest	41.5"	44"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES	YES
Warranty on Electrical Parts	3 YR	3 YR
Pro Rated Years 4 Thru 7	YES	YES
Warranty on Mechanical Labor	3 YR	3 YR
Warranty on Electrical Labor	1 YR	1 YR
Hand Control Routing	Inside Pocket	Inside Pocket
Extra Pocket Option	Standard	Standard
Heat & Vibrating Massage Option	NO	NO
HeatWave Technology Option	NO	NO
Nirvana with HeatWave Option	NO	NO
Articulating Headrest Option	NO	NO
Power Adjustable Lumbar Option	NO	NO
Foot Rest Extension Option	NO	NO
Straight Lift Option	Standard	Standard
USB Charging Port	Hand Control	Hand Control
Battery Backup	9 volt	9 volt
SmartTek Diagnostic System	Standard	Standard
Place of Origin	Imported	Imported



Shown in Calypso


Comfort Zone™ by GOLDEN

Comfort Zone 1 provides ergonomic placement of your lower body allowing your muscles to relax for long-lasting comfort. Comfort Zone 1 can be adjusted using the two buttons on the hand control.



LOWER BODY UPPER BODY HEADREST LUMBAR TWILIGHT TILT

WE PATENTED COMFORT™



Scan with your smartphone camera to learn more about Comfort Zone by Golden.



Phone: 800-624-6374
Fax: 800-628-5165
www.goldentech.com

FDA Class II Medical Devices designed by Golden Technologies are intended to assist persons with physical disabilities.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Complies with
UL 73 and
CSA C22.2 No. 68



www.facebook.com/goldentechnologies

LIT-SHT-PR401, Rev. B 010925