



**Rhea
PR442**



Comfort Zone™
by GOLDEN

Take a vacation from your busy day in the Rhea. Its overstuffed seam back, plush chaise pad, and open arm design offers a cup holder, storage compartments, power port with removable LED light and tray table, wireless phone charger and a USB port for convenience. Enjoy customizable comfort with our 4-Zone Comfort Zone system that includes our patented MaxiComfort® Positioning. For extra wellness and comfort, our soothing Far-infrared warmth awaits with HeatWave™ Technology, delivering therapeutic heat from shoulders to feet for ultimate relaxation.




FDA Class II Medical Devices designed by Golden Technologies are intended to assist persons with chronic physical disabilities.

Standard Fabrics:



Lagoon



Umber

The Sinda polyurethane pattern is soft and supple and was designed to be breathable yet highly durable. It provides a realistic leather-like look with a subtle texture. The composite material and polyurethane resin coating make this fabric water resistant, flexible, and lightweight.

Crossover fabrics are unavailable.

Limited Edition Fabric:



Cork*

100% Top Grain Italian Leather!

*Upcharge applies.

Standard Features:

- Patented MaxiComfort® Positioning
- Adjustable Headrest & Lumbar
- Programmable AutoDrive™ Hand Control
- HeatWave™ Technology
- Wireless Smartphone Charger
- USB Charging Port
- Built-in Cup Holder
- Open Arm Storage
- SmartTek® Diagnostic System
- Extra Pocket
- Power Port with Removable LED Reading Light & Tray Table

**NEW Power Port
with Optional
Accessories!**

DeLUNA Series PR442

Model Number	PR442-MED	PR442-LAR
Size	Medium	Large
Class II Medical Device	YES	YES
Weight Capacity	400 LBS	400 LBS
Positioning	MaxiComfort	MaxiComfort
Comfort Zone(s)	4	4
Twilight Technology	NO	NO
ZG+	NO	NO
Chaise Seat	YES	YES
Pocketed Coil Seat	YES	YES
Back Type	Seam	Seam
Removable Back Type	Quick Release	Quick Release
Head Rest Cover	NO	NO
Arm Rest Covers	NO	NO
Overall Width	37"	39"
Overall Height	45"	50"
Seat Width	22"	24"
Floor to Top of Seat	20"	21"
Seat Depth	21"	22"
Seat to Top of Back	29"	33"
Seat to Top of Arm	7"	8"
Distance Required from Wall Reclined	25"	30"
Extension Reclined from Back of Seat to End of Foot Rest	43"	44"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES	YES
Warranty on Electrical Parts	3 YR	3 YR
Pro Rated Years 4 Thru 7	YES	YES
Warranty on Mechanical Labor	3 YR	3 YR
Warranty on Electrical Labor	1 YR	1 YR
Hand Control Routing	Right Arm Compartment	Right Arm Compartment
Extra Pocket Option	Standard	Standard
Heat & Vibrating Massage Option	NO	NO
HeatWave Technology Option	Standard	Standard
Nirvana™ with Heatwave Option	NO	NO
Articulating Headrest Option	Standard	Standard
Power Adjustable Lumbar Option	Standard	Standard
Foot Rest Extension Option	NO	NO
Straight Lift Option	Standard	Standard
USB Charging Port	Left Arm Compartment	Left Arm Compartment
Battery Backup	Lithium-Ion	Lithium-Ion
SmartTek Diagnostic System	Standard	Standard

Designed and engineered in Old Forge, PA. Assembled in China using Golden components.

Standard Features & Accessories:



Cup Holder, Power Port, USB Charging Port & Storage Compartment



Wireless Smartphone Charger & Storage Compartment



Adjustable Headrest & Lumbar



HeatWave™ Technology



Removable LED Reading Light & Tray Table

Optional Accessories:

- Cooling Fan
- Flex Light
- Audio Torch
- Illuminating Snack Bowl
- Essential Oil Diffuser
- Tablet Holder
- Stem Glass Holder




Comfort Zone™
by GOLDEN

Comfort Zones 1 - 4 provide an ergonomic comfort experience with an emphasis on head and back support. Each Comfort Zone can be adjusted independently by using the exclusive AutoDrive hand control.



WE PATENTED COMFORT™



Scan with your smartphone camera to learn more about Comfort Zone by Golden.



Phone: 800-624-6374

Fax: 800-628-5165

www.goldentech.com

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Complies with
UL 73 and
CSA C22.2 No. 68



www.facebook.com/goldentechnologies

LIT-SHT-PR442, Rev. A 070125