




## Standard Fabrics:



Antler




Elk



Lakefront

This contemporary plush fabric comes in beautiful, rich colors. It features our exclusive SofTouch™ with Stain Defense coating.



Chestnut

Our Urethane fabric is highly durable and provides the most realistic imitation of leather. The composite material and polyurethane coating make this fabric water resistant, flexible and lightweight.

Crossover fabrics are unavailable.

**Elara**  
**PR118**



*Comfort Zone™*  
by GOLDEN

Inspired by the moon, the DeLuna™ Series of luxury power lift recliners is designed to offer you the most calming comfort and restorative sleep. The Elara is offered in two sizes with a back style that provides maximum comfort in a tasteful seam design. It features a supportive chaise pad and a USB port on the hand control. Available in luxurious new fabrics to match any décor. Enjoy the powerful energy and take time to better connect with yourself in your DeLuna recliner by Golden!



## Standard Features:

- Easy Grip Armrest
- Soft-Feel Hand Control
- USB Charging Port



## DeLUNA Series PR118

Model Number	PR118-MSM	PR118-LAR
Size	Med-Small	Large
Weight Capacity	400 LBS	400 LBS
Positioning	3	3
Comfort Zone(s)	1	1
Twilight Technology	NO	NO
Chaise Seat	YES	YES
Pocketed Coil Seat	NO	YES
Back Type	Seam	Seam
Removable Back Type	Bracket	Bracket
Head Rest Cover	YES	YES
Arm Rest Covers	NO	NO
Overall Width	32"	34"
Overall Height	43"	46"
Seat Width	21"	23"
Floor to Top of Seat	18.5"	19.5"
Seat Depth	21"	21"
Seat to Top of Back	29"	31"
Seat to Top of Arm	7"	6"
Distance Required from Wall Reclined	19"	20"
Extension Reclined from Back of Seat to End of Foot Rest	42"	42"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES	YES
Warranty on Electrical Parts	3 YR	3 YR
Pro Rated Years 4 Thru 7	YES	YES
Warranty on Mechanical Labor	3 YR	3 YR
Warranty on Electrical Labor	1 YR	1 YR
Hand Control Routing	Inside Arm	Inside Arm
Extra Pocket Option	NO	NO
Heat & Vibrating Massage Option	NO	NO
HeatWave Technology Option	NO	NO
Nirvana with HeatWave Option	NO	NO
Articulating Headrest Option	NO	NO
Adjustable Lumbar Option	NO	NO
Foot Rest Extension Option	NO	NO
Straight Lift Option	YES	YES
USB Charging Port	Hand Control	Hand Control
Battery Backup	YES	YES
Place of Origin	China	China



Elegant seamed back design and supple urethane fabric make this chair both pleasing to the eye and highly durable.



Shown in Chestnut

Designed and engineered in Old Forge, PA. Assembled in China using Golden components.

*Comfort Zone*  
by GOLDEN

**Comfort Zone 1** provides ergonomic placement of your lower body allowing your muscles to relax for long-lasting comfort. Comfort Zone 1 can be adjusted using the two buttons on the hand control.



**WE PATENTED COMFORT™**


FDA Class II Medical Devices designed by Golden Technologies are intended to assist persons with physical disabilities.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Complies with  
UL 73 and  
CSA C22.2 No. 68



Scan with your smartphone camera to learn more about Comfort Zone by Golden.



Phone: 800-624-6374  
Fax: 800-628-5165  
[www.goldentech.com](http://www.goldentech.com)



[www.facebook.com/goldentechnologies](http://www.facebook.com/goldentechnologies)

LIT-SHT-PR118, Rev. B 092625