

Standard Fabrics:



Acorn



Graphite



Maple

This faux leather is a rich and modern alternative to natural leather that is extremely comfortable while offering the rugged durability of leather. Made from 100% polyester, the fabric is easy to care for ensuring functionality as well as elegance. This long-lasting fabric will complement a variety of decors while standing the test of time under heavy use.

Crossover fabrics are unavailable.

Dione PR446





The Dione is offered in both medium and large sizes and features a stylish, overstuffed seamed back and plush chaise pad. It comes standard with adjustable lumbar and headrest, programmable AutoDrive™ hand control, and USB charging port for your phone or tablet. Also standard is our patented four Comfort Zone positioning system, which provides a stress-free and supportive reclining experience.

Standard Features:

- Patented MaxiComfort® Positioning
- Adjustable Headrest
- Adjustable Lumbar
- Programmable
 AutoDrive™ Hand
 Control
- USB Charging Port
- Extra Pocket
- SmartTek® Diagnostic System





DeLUNA Series PR446

| Model Number | PR446-MED | PR446-LAR |
|---|-------------|---------------|
| Size | Medium | Large |
| Class II Medical Device | YES | YES |
| Weight Capacity | 400 LBS | 400 LBS |
| Positioning | MaxiComfort | |
| Comfort Zone(s) | 4 | 4 |
| Twilight Technology | NO | NO |
| ZG+ | NO | NO |
| Chaise Seat | YES | YES |
| Pocketed Coil Seat | YES | YES |
| Back Type | Seam | Seam |
| Removable Back Type | | Quick Release |
| Head Rest Cover | NO | NO |
| Arm Rest Covers | NO | NO |
| Overall Width | 34" | 35" |
| Overall Height | 43.5" | 46" |
| Seat Width | 22" | 22.5" |
| Floor to Top of Seat | 20.5" | 22" |
| Seat Depth | 21" | 21.5" |
| Seat to Top of Back | 28" | 30" |
| Seat to Top of Arm | 7" | 7.5" |
| Distance Required from Wall Reclined | 24" | 26" |
| Extension Reclined from Back of Seat | | |
| to End of Foot Rest | 42.5" | 44" |
| Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism | YES | YES |
| Warranty on Electrical Parts | 3 YR | 3 YR |
| Pro Rated Years 4 Thru 7 | YES | YES |
| Warranty on Mechanical Labor | 3 YR | 3 YR |
| Warranty on Electrical Labor | 1 YR | 1 YR |
| Hand Control Routing | Inside Arm | Inside Arm |
| Extra Pocket Option | Standard | Standard |
| Heat & Vibrating Massage Option | NO | NO |
| HeatWave Technology Option | NO | NO |
| Nirvana with HeatWave Option | NO | NO |
| Articulating Headrest Option | Standard | Standard |
| Power Adjustable Lumbar Option | Standard | Standard |
| Foot Rest Extension Option | NO | NO |
| Straight Lift Option | Standard | Standard |
| USB Charging Port | Right Arm | Right Arm |
| Battery Backup | Lithium-ion | Lithium-ion |
| SmartTek Diagnostic System | Standard | Standard |
| Place of Origin | Imported | Imported |





Designed and engineered in Old Forge, PA. Assembled in China using Golden components.





Comfort Zones 1 - 4 provide an ergonomic comfort experience with an emphasis on head and back support. Each Comfort Zone can be adjusted independently by using the exclusive AutoDrive hand control.





UPPER BODY









WE PATENTED COMFORT™





Sizes may vary depending on fabric, filling material, upholstery and carpet

thickness. All measurements taken on concrete floor using levels and metal rulers.





Phone: 800-624-6374

Fax: 800-628-5165 www.goldentech.com

Scan with your smartphone

camera to learn

more about

Comfort Zone by Golden.